**DESTINATION DU DISPOSITIF**

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
## User manual HR foam APLOT alveolar mattress

### DEVICE INTENDED USE

This high resilience foam static mattress is intended to be used in combination with a medical bed for specific pressure ulcer prevention purposes due to a loss of autonomy and/or deterioration of health resulting in reduced mobility and the risk of pressure ulcers.

### INDICATIONS

Bedridden patients with reduced mobility and "Low to Medium" risk of pressure ulcers (according to the Braden scale or any other validated scale and on medical advice). Assists in the treatment of stage I to 2 pressure sores according to EPUAP, up to stage 4 outside the support area.

### CONTRAINDICATIONS

Patient weight level below or above the limits defined in the technical data table. For use on medical imaging table.

### ADVERSE SIDE EFFECTS

Not documented.

**Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.**

### MEDICAL DEVICE COMPOSITION/TECHNICAL DATA

References Models and versions*	VA103MB14	VA103MB14RC	VA103MB14/HIP	VA103MB14/HIPH
Sizes	195x85x14 cm		195x85x14 cm	
Lower specs	HR: 42.9 kg/m <sup>3</sup> 4.4 kPa		HR: 42.9 kg/m <sup>3</sup> 4.4 kPa	
Patient weight min and max	40 to 120 kg		40 to 120 kg	
Foam material	Dermalon FR (Fire resistant)		100% coated polyurethane on PES weave	

\*The indicated references may be completed with a suffix to constitute a commercial reference depending on the configuration (type of protective cover, packaging) and/or medical equipment dealer

### CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION

**DEVICE PERFORMANCE CHARACTERISTICS:** By immersion, this support relieves pressure on the body by improving blood circulation and thus reducing the risk of pressure ulcers on the supporting areas of the body.

**EXPECTED CLINICAL BENEFITS:** Maintenance of tissue oxygenation at the cutaneous and subcutaneous tissue level in contact with the support.

**INFORMATION FOR HEALTHCARE PROFESSIONALS:** Monitor the condition of the patient's skin in contact with the support several times a day. Use pressure relief devices or positioning systems in patients with pressure ulcer(s).

### PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE

**TRAINING AND QUALIFICATION OF THE DEVICE USER:** user training must be carried out by individuals trained and approved by the economic operators concerned, notably in terms of safety and non-compliance reporting.

**DEVICE installation:** The mattress is delivered with its draw sheet or removable zipped cover. The product is ready to be installed. The surface of the foam composed of studs ("alveolar") must be positioned in contact with the patient's body. The protection must be dry.

**Preventive maintenance:** Carry out a regular visual inspection of the foam condition: visible sagging of material and a non-consistent and very slow return of the foam shape are ageing criteria that compromise the properties of the support. Annual inspection of the condition of the cover (surface aspect and zippers) by exposing the inner surface to a light source: no holes and/or tears. Replace the cover if its surface appearance changes.

### CLEANING AND DISINFECTION

Promust cover PU and HDPU						
	Moderate washing up to 90°C	Maximum chlorine concentration of 500ppm	Do not iron	Do not dry clean	Dry with reduced thermal restriction	Use an approved surface disinfectant/detergent
Dermalon draw sheet						

### INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE

The product must exhibit physical and bacteriological cleanliness. "Promust" protection treated with a biocide substance at no risk to users. Avoid wet contact with the foam. Prohibit all scouring, stripping products or solvents and "piercing-sharp" objects in direct contact with the cover.

### WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED

#### PRECAUTIONS FOR USE

A bedsore is a lesion of varying depth in the skin, associated with excessive and prolonged compression of tissue between the body and a support surface. This excessive pressure can prevent blood circulation and lead to pressure ulcers, a dangerous lesion because it develops from the inside of the body to the outside, and is therefore initially invisible. A pressure sore can have several forms: a simple persistent rash for more than a day, hardening of the skin, a wound of varying depth that can in severe cases affect muscles or underlying bone. Pressure (max) in mmHg before ageing, at the maximum validated weight, measured flat (0°)

SUITABILITY FOR USE	Max weight User	Head	Buttocks	Heel
APLOT mattress (mm of hg)	140 kg	54	46	92

#### WARNING

(Re)assess the risk of patients being trapped in the non-moving parts of the medical bed associated with the "therapeutic mattress and accessories and the articulated bed base positions". The use of segmented or full length side barriers must be required.

External cardiac massage not compatible without a board placed between the thorax and the top surface of the mattress

#### ACTION REQUIRED

- A mattress must be placed on a good quality, flat, spring mattress. An over-mattress must be used on a mattress in good condition.
  - Make sure the mattress size is compatible with the bed. The system must be used with its original protection.
  - When applicable, check that the bed railing height is compatible with the mattress thickness or mattress and over-mattress thickness (at least 22 cm between the top of the mattress and the top of the railings)
  - Important: a support itself is not enough to prevent pressure ulcers; other preventive measures are also essential:
    - change position frequently (at least every 2 to 3 hours);
    - maintain skin hygiene and avoid maceration; do not tuck in bedding tightly and use elastic sheets to keep mattress effectiveness.
    - in the event of incontinence, changing protection regularly;
    - examine skin condition (or have it examined) every day;
    - make sure food intake is sufficient and drink regularly and in sufficient quantities.
- If any of these measures cannot be followed, your doctor or nurse must be informed as quickly as possible.
- For the support to be effective, it is important to limit over-thickness between the body and the support as much as possible, except for the sheet for a bed support, clothing and a potential complete change. Favour loose cotton clothing and if possible without stitching in the support zone. Do not insert: folded towels or sheets, additional cushions, etc.
  - Ensure there are no foreign objects such as: tubes, crumbs, fat, etc.
  - Important: consult your prescriber within one month of purchasing this support.

#### CIRCUMSTANCES IN WHICH USERS SHOULD CONSULT WITH A HEALTHCARE PROFESSIONAL

Immediately notify your doctor or nurse of all abnormal events such as fever, pain or redness or whitening of pressure points (head, shoulder, back, hip, shoulder blade, pelvis, heel, etc.).

#### INFORMATION ON ALL KNOWN RESTRICTIONS ON THE COMBINATION WITH SYSTEMS AND EQUIPMENT

The use of accessories other than those provided and specified by WINNOCARE may result in the system malfunctioning and/or render it incompatible with safe use and the preservation of its performances.

#### STORAGE, HANDLING, DISPOSAL

##### STORAGE AND USAGE CONDITIONS

The supports must ideally be stored flat, away from direct light and excessive humidity.

	Use	Storage
Temperature range	+15°C to +45°C / +59°F to +113°F	-25°C to 70°C / +13°F to +158°F
Humidity range	30% - 70%	30% - 95%
Atmospheric pressure	50kPa - 106 kPa // 7.25 Psi / 15.37 Psi	

##### ROLLED MATTRESSES

APLOT: 30 mn → 2h : This time may be increased if the storage temperature is lower than 15°C and/or in line with the storage period (0 to 5 months).

The return to a normal visual appearance of the mattress under normal usage temperature conditions may take up to 24 hours for ALOVA, APLOT and MORFÉA mattresses but this does not present any risk for their use. Small creases may be visible on the lower surface of the CIC cover after unrolling. They present no risk to the patient and do not affect the efficiency of the final product.

##### SERVICE LIFE

Highly resilient mattresses have an estimated service life of 6 years.

The rolled-compressed mattress life span is 6 months (see date on package label).

##### PRODUCT DISPOSAL

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

##### GUARANTEE

APLOT mattresses have a 2-year guarantee. This guarantee starts from the date of purchase of the product from your distributor. This guarantee does not cover normal wear and tear of the product and its protection. Please contact your distributor about a defective product. They will carry out the necessary procedures with our company to proceed with either repair or exchange.

##### Class I medical device according to Regulation (EU) 2017/745.

Management systems: Quality: ISO 13485:2016 - Environmental: ISO 14001:2015

Technical compliance certificates (FCBA): ESC N° 16-005 of 10 June 2016

#### Contents of package

- 1 mattress
- 1 protective cover
- 1 user manual

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