



Archimedes Bath Lift

The Archimedes Bath Lift will lower a person into a bath, on a fixed position seat, using a battery powered hand control.

Easy to fit

- Simple to fit – separates into three pieces, making it portable and lightweight
- Fits further back in the bath to maximise leg room and to avoid most bath handles and anti-slip coatings
- Unique 4 way quick release suckers make removing the Bath Lift simple – reducing strain on your back

Simple to maintain

- The easy to clean steel base unit has been designed with a structure that maintains strength whilst reducing weight
- No LOLER test required for two years
- Divides into three so it can be removed and cleaned with ease

Safe to use

- Unique flexible flaps that bend on contact with an obstruction, limiting the damage caused by entrapment
- Bath flaps have an anti-pinch design, eliminating the risk of skin being trapped between the seat and the flap on transfer
- Only lowers a person into the bath when there is sufficient battery power to raise them back out again

Improves well-being



- Bathing can reduce pain and promote a feeling of well-being
- Equipment will reduce the amount of moving and handling required from a carer
- Can be removed from the bath so enabling a more flexible environment for shared bathrooms

Technical Specifications

Max user weight	165 kg
Max height	45 cm
Min height	8 cm
Seat depth	44.5 cm
Total length	61 cm
Width (ex side flaps)	39.5 cm
Width of side flap	17 cm
Total product weight	10 kg
Back/headrest length	66.5 cm



Scan me



Product code: LAA3702A

Features

1. **Two anti-pinch transfer flaps**
2. **Sealed Hand Control** with a handset which has easily visible control buttons.
3. **Release the transit tie** by pulling up the 'O' ring and lifting the ball clear of the slot in the locating collar on the actuator.
4. **Battery powered by actuator**
5. **Strong suckers** which keeps the Archimedes in place whilst bathing.
6. **Safety feature** which ensures there is enough charge to lift user out of bath.

Accessories

- Archimedes Spare Pack
- Archimedes Hand Control
- Archimedes Charger
- Archimedes Side Flap Protector Set

