

# MÉDI REPOS® - MÉDI DETENTE®



MÉDI REPOS CHAIR with footstool.  
Banana Yellow PVC, Structure structure ral 9016

## MÉDI REPOS

## WOODEN CHAIR

tilting up to 40°



### CONVENIENT

- ▶ Back can tilt up to 40° using a gas spring.
- ▶ Handles for tilting the back on both sides.
- ▶ Small mobility castors on the 2 rear legs. (movement when empty)

### STYLISH

- ▶ The wooden parts, in solid beech, are stained or lacquered to protect the wood and make it easier to clean.
- ▶ Choice of colors according to our color chart (cf. color chart page).
- ▶ PVC coating according to our color chart (cf. color chart page).

## MÉDI DETENTE

## METAL CHAIRS

### CONVENIENT

- ▶ Back can be tilted by gas spring up to 40°.
- ▶ Back 2/3 1/3.
- ▶ Push bar.
- ▶ Tilting handles on the side.
- ▶ Retractable armrests, adjustable into 4 positions.
- ▶ Removable non-slip seat cushion.
- ▶ Chair covers can be separated and removed.



FOOTSTOOL  
Leg rest.  
Round tube

### SAFETY

- ▶ 4 x 100 mm castors including 2 with brakes.
- ▶ Thermoformed footrest plate that can be lifted with anti-tip mechanism.

### ACCESSORIES

- ▶ IV stand, tray...
- ▶ Epoxy structure, RAL 7035 grey (one color)
- ▶ PVC coating according to our color chart (cf. color chart page).



## MÉDI OCÉAN

- ▶ Back can be tilted by gas spring (40°).
- ▶ Retractable armrests.
- ▶ Round tube.
- ▶ Tilting handles on both sides.
- ▶ Height adjustable seat.

## ERGONOMICS

- ▶ MÉDI REPOS chairs facilitate the transfer of patients from bed to chair and allows them to rest comfortably.

## ROBUST

- ▶ Structure in solid wood or epoxy metal.

## TECHNICAL SPECIFICATIONS

DESCRIPTION	SEAT				BACK Height	ARMRESTS Height	OVERALL		
	Width (between armrests)	Depth	Seat height from ground	Cushion thickness			Width	Height	Depth
WOODEN CHAIR									
MÉDI REPOS	52	48	46	10	73	66/59	69	117	74
FOOTSTOOL	43	47	45	10	NA	NA	48	45	51
METAL CHAIRS									
MÉDI OCÉAN	55	48	46.5 / 50.5 / 54.5	10	73	66	72	117	79
MÉDI DETENTE	49	45	50	10	73	50 / 69	67	120	79
FOOTSTOOL	43	47	45	10	NA	NA	50	45	48

Patient weight supported: 135 kg

warranty: 2 years

Photos not contractually binding

FP0000026/03-2017

