

Mobile Flex

4-post system

Full adjustable, room covering, free-standing traverse system

Primary application

Home care
Institutional
Hospital

Room covering and free standing traverse system from Ergolet.

Mobile Flex is a flexible system that is easy to install and can be used right after the installation. Mobile Flex is installed and dismantled without any changes to the room.

Adjusting the height via the gas springs is carried out in just a few seconds

Due to its many adjustment options, Mobile Flex meets all lifting and transferring needs.

Combined with the Luna ceiling lift, Mobile Flex is a simple and functional lifting system.

Key standard functions

Max load of 275 kg
Applies full capacity of the Luna ceiling lift
Available in several variants
Covers up to 4x4 meters
Fully adjustable 4-pole system
Bayonet trolley
Can be used outdoors

Accessories

Multiswitch



PRODUCT HIGHLIGHTS

- Practical design and simple functionality
- Room covering system
- Easy to clean



Easy to adjust the width without tools

Version	Item-number
Mobile Flex 4-post system: model 3x3 meter	016-01220
Mobile Flex 4-post system: model 3x4 meter	016-01221
Mobile Flex 4-post system: model 4x3 meter	016-01222
Mobile Flex 4-post system: model 4x4 meter	016-01223



Easy height adjustment (gas spring mechanism). Simple user manual is located above the handle.


TECHNICAL SPECIFICATION

Max. load 275 kg

	3x3 (track x traverse)	3x4 (track x traverse)	4x3 (track x traverse)	4x4 (track x traverse)
a	1891 mm	1891 mm	1891 mm	1891 mm
b	2160 mm	2160 mm	3160 mm	3160 mm
c	3160 mm	3160 mm	4160 mm	4160 mm
d	2550 mm	2550 mm	2550 mm	2550 mm
e	2330 mm	3330 mm	2330 mm	3330 mm
f	3330 mm	4330 mm	3330 mm	4330 mm



Luna ceiling hoist mounted on the Mobile Flex



Consult the Technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Ergolet reserves the right to changes without notice.

Product photos may include accessories. Consult your dealer.

CE MDD93/42/EEC
Class I

WinnCare Nordic ApS
Taarnborgvej 12C
4220 Korsoer
Denmark
tel. +45 70 27 37 20
e-mail: info@winnCare.dk

WINCARE
Nordic

date 09.10.2017 - ver. 7

www.winnCare.dk