

# POOL E-185<sup>®</sup>




## FOR PATIENTS WEIGHING UP TO 185 KG

### SAFETY

- ▶ Control box with emergency stop and built-in load indicator.
- ▶ Audible alarm: to warn when battery is running low.
- ▶ Electric emergency lowering system.
- ▶ Nylon coated stainless steel metal structure. Control box with emergency stop and built-in load indicator.
- ▶ Audible alarm: to warn when battery is running low.
- ▶ Electric emergency lowering system.
- ▶ Nylon coated stainless steel metal structure.

### CONVENIENT

- ▶ Designed for use in swimming pools.
- ▶ 2 rotation bases:
  - To be cast in concrete
  - To be screwed directly into the ground
- ▶ Foldable
- ▶ 360° spreader bar rotation
- ▶ Accessories:
  - Special stretcher for swimming pool use
  - Special chair for swimming pool use



Ref: 102195001:  
rotation base to be cast  
in concrete.




Ref: 102195002:  
rotation base to be screwed  
directly into the ground.



Ref: 102106004



Ref: 102106003







Detachable battery

### THE BATTERY

can be detached and charged up on the hoist or with an external charger (optional).

- 1 Charge indicator.
- 2 Stop push button.
- 3 Emergency electric lowering.



CONFORMS TO STANDARD EN ISO 10535: 2007 APPLICABLE TO HOISTS FOR TRANSFERRING PEOPLE WHO ARE DISABLED.

### ADAPTED HOIST

- ▶ Hoist is suitable for patients weighing up to 185 kg.

### TRANSFER

- ▶ Transfer from wheelchair to swimming pool.

### ERGONOMIC

- ▶ Lift from the ground.
- ▶ elevation range up to 2.10 m.
- ▶ Wide handles.

### ROBUST

- ▶ LINAK motorization and materials.

DESCRIPTION	REFERENCE	MAX PATIENT WEIGHT	STANDARDS	MOTORIZATION
POOL E-185	102102006	185 kg	ISO 10535 EN 12182	LINAK IP 66

Remote control, electric and mechanical parts: 5-year warranty  
Batteries, slings and weighing unit: 1 year