



e-STATION IN-TRACK CHARGING SOLUTION

e-Station, the newest addition to the family of WinnCare track solutions, ensures the Luna ceiling hoist is always charged and ready for duty. Simply just park Luna in the charging area of the e-Station, and charging will automatically begin.



An aesthetic solution

The e-Station is a 1 meter piece of 75 mm track with an integrated charging feature.

The Luna hoist is aesthetically pleasing and our new e-Station compliments by offering a discreet solution free of a visible charging station mounted on the wall.

For both new installations and retrofitting

The e-Station can easily be mounted and implemented in new installations and also retrofitted into existing single track systems.

It's compatible with turntables and with the right configuration of tracks, it can also be fitted into XY-systems.

The e-Station can be either end mounted on wall or suspended from ceiling.

Two Lunas at a time

e-Station is compatible with Luna 200/275 kg., fitted on a bayonet cart with power transfer.

Integrating the e-Station in an existing single track system is simple and requires only a change of the Luna side foil and cart; two easy steps and the Luna is ready for charging in the e-Station.

The e-Station's charging zone fits two Lunas at a time.