


ALOVA™ XL



Braden scale score



KEY FEATURES

- ▶ Single-piece mattress made up of:
 - 7 cm molded visco elastic top layer: density: $\geq 75 \text{ Kg/m}^3$; hardness: 1.7 kPa
 - 10 cm polyether foam base: density: 28 Kg/m^3 ; hardness: 3.6 kPa
- ▶ Patient weight: 30 to 150 Kg
- ▶ Dimensions: VMA117XL: 196 x 117 x 17 cm (77.1 x 46 x 6.6 inches)
VMA137XL: 196 x 137 x 17 cm (77.1 x 53.9 x 6.6 inches)
- ▶ Mattress weight: 24 Kg (VMA117XL) and 28 Kg (VMA137XL)
- ▶ Cover available: Promust PU

KEY BENEFITS

- ▶ Visco elastic foam provides excellent pressure re-distribution properties: homogeneous immersion of the body into the foam, thus increasing the bearing surface and reducing pressure applied in areas at risk of pressure ulcers.
- ▶ Suitable for large width sleeping surfaces of 120 or 140 cm width medical beds.
- ▶ Visco elastic foam provides patient comfort.

INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT MEDIUM TO HIGH RISK OF PRESSURE ULCER DEVELOPMENT, BEDRIDDEN MORE THAN 15 HOURS A DAY.

HELP IN TREATING EXISTING PRESSURE ULCERS FROM STADE 1 TO 2 COMBINED WITH TECHNICAL AID DEVICES.

PATIENT WEIGHT: 30 TO 150 KG/ 4.7 TO 23.6 ST

- ▶ Single-piece mattress made up of visco elastic foam top layer on High Resilient (HR) polyurethane foam
- ▶ ALOVA XL: therapy and comfort in large width bed
- ▶ Removable and washable integral cover

DESCRIPTION	ITEM REFERENCES ACCORDING TO PACKAGING		FIRE STANDARDS
	CARDBOX	ROLLED AND COMPRESSED (RC)	
ALOVA™ XL 120 cm width, Promust PU cover	VMA117XL	VMA117XL/RC	EN 597-1 and 2
ALOVA™ XL 140 cm width, Promust PU cover	VMA137XL	VMA137XL/RC	

3 years warranty

Cleaning instructions: Promust PU cover



FP000037/01.2018

Photos not contractually binding

